

Opp Scan Help

Printer Commands

- **Connect**
Attempts to connect to the specified baud rate and port.
- **Add/Subtract X (+/-)**
Adds or subtracts the value specified in the spinbox to the current value of X.
- **Add/Subtract Y (+/-)**
Adds or subtracts the value specified in the spinbox to the current value of Y.
- **Add/Subtract Z**
Adds or subtracts the value specified in the spinbox to the current value of Z.
- **Baud Rate**
The measure of the speed of data transmission between two devices.
 - **Mini printer:** 9600
 - **I3:** 115200
- **Conductance Mode**
Sets the sample mode to use conductive touch.
 - **Note:** You must have started conduction to use this mode.
- **Dimension Mode**
Circle beside the dimension in the sampling.
 - Disables the dimension and calculates it using resolution and sampling spots.
- **Disconnect**
Disconnect from the 3D printer.
- **Dwell Time**
The time in between sampling spots where it sits.
- **Pause**
Pauses sampling and sends the node to X10 Y10 Z100.
- **Port**
Place within an operating system network where connections start and end.
 - **Example:** The mini printer will be called "Original Prusa MINI" on Mac and will appear as a COM port on Windows.
- **Resolution**
Distance between each sampling location in either the X or Y direction.
- **Resolution Mode**
Disables the resolution and calculates it using the dimension and sampling spots.
- **Sample Time**
The amount of time in seconds that the probe will pause after it hits the lower Z or detects conduction.
- **Sample Z**
The position in the Z axis that the program will sample during constant Z mode. Used for flat surfaces.

- **Example:** For dewaxing, use a sample Z equal to your start Z for a continuous stream.
- **Set Z (Sample Z)**
Sets the sample Z to the current Z position.
- **Sampling Spots**
The number of spots the sampling program will check in the respective X and Y directions.
- **Sampling Spots Mode**
Disables sampling spots and calculates using resolution and dimensions, rounding to the nearest significant number (e.g., 3.2 is rounded to 4).
- **Set All**
Sets the start position to the current position for each coordinate.
- **Set X, Set Y, Set Z**
Sets the start position to the current position for the X, Y, or Z coordinates, respectively.
- **Start Position**
The start position for each coordinate when sampling begins. The program will return to this position between each sampling spot.
- **Start Run**
Starts the sampling program with the current G-code.
- **Step Size**
The size of each command sent in conductance mode in the Z direction.
 - **Note:** Conductance only triggers a stop after a command has been completed. A larger step size is generally faster.
- **Stop**
Stops the sampling and sends the printer to X10 Y10 Z100.
- **Update**
Updates all the ports to check for any new connections.
- **XY Speed**
Speed at which the sampling moves in the XY position in mm/min.
- **Z Down Speed**
The speed at which the probe moves down during sampling for both conductance and constant Z (in mm/min).
- **Z Up Speed**
The speed at which the probe moves up (in mm/min).

Temperature Controls

- **Set Temperature**
Sets the temperature of the extruder and the bed as specified by the respective spinboxes.
- **Bed Temperature**
Temperature of the bed in Celsius.
- **Extruder Temperature**
Temperature of the extruder hot-end in Celsius.

Conductance Commands

- **Connect**
Connects to the specified baud rate and port.
- **Disconnect**
Disconnects the conductance.
- **Baud Rate**
The measure of the speed of data transmission between two devices.
 - **Arduino:** Usually 9600.
- **Port**
Place within an operating system network where connections start and end.
 - **Mac:** Port will start with a serial connection name.
 - **Windows:** Will be a COM port.
- **Update**
Updates all the ports to check for any new connections.
- **Zero At**
Zeros the conductance graph at a specific value.
 - **Example:** Setting the value to 100.
 - 99 -> 0
 - 100 -> 0
 - 101 -> 1
- **Threshold**
Used in conductance mode to specify the minimum conductance to stop at.
- **Start**
Starts the conduction graph and resets it if needed.
- **Stop**
Stops the conduction graph.
- **Save**
Prompts the user to enter a file to save the conductance data under.
 - **Saved columns order:** Time (ms), conductance, X coordinate, Y coordinate, Z coordinate.
- **Scale**
Adjusts the displayed size of the conductance graph.

Camera Controls

- **Connect**
Connects to the specified camera.
- **Disconnect**
Disconnects from the camera.
- **Capture**
Starts the camera.
- **Camera**
Select the camera to connect to.

- **Refresh Rate**
The number of times the image on the camera's display is refreshed per second.
- **Take Photo**
Takes a photo and labels it with the current time.
- **Start/Stop Video**
Starts and stops video recording.
 - **Saves as** `<time>.mp4`.
- **Pause**
Pauses the video recording.

Pump Controls

- **Baud Rate**
The measure of the speed of data transmission between two devices.
- **Connect**
Attempts to connect to the pump under the specified port and baud rate.
- **Constant Flow Rate**
Select to run with a consistent flow rate.
- **Disconnect**
Disconnects from the pump.
- **Dwell Time**
Time at dwell position.
- **Flow Rate**
The quantity of fluid passing through a cross-section of a pipe in a specific period of time.
- **Flow Rate Program**
Additional flow rate control options.
- **Port**
The port used for the pump connection.
- **Sampling Rate**
The rate at which samples are taken.
- **Set Flow Rate**
Sets the flow rate for the pump.
- **Set Size**
Sets the size for the pump operation.
- **Set Volume**
Sets the volume for the pump operation.
- **Start**
Starts the pump operation.
- **Stop**
Stops the pump operation.
- **Syringe Size**
Specifies the size of the syringe used.

- **Volume**

Specifies the volume for the pump operation.

-